

Customizing the CDEMS progress note

The user can define any condition and related diagnoses, meds, services and lab in CDEMS.

Here is the procedure for customizing CDEMS:

Access the Setup wizard from the opening screen of the CDEMS data entry program by pressing **F11** function key on your keyboard.

It will take you to database view> click on **Forms**> double click on **frmSetup**.

It will open new window, which is yellow green in color.

In the green box where it says Step 1 select DM-2 or DM-1 for health condition.

Note: If you want to add a new health condition just type it in the drop down menu and you will be prompted to add the new condition.

Remember to always select a health condition before using the setup wizard.

In the box labeled New Rx or New Services or New Lab type the appropriate new med or new service or new lab. Press enter. Follow the instructions on the screen.

The new variable will be added in the right- sided box of database of all stored values.

Highlight the newly created variable and Use the “< “ to move the new variable to the left box.

The contents of the left box are the one that you see on your progress notes.

If you want to take out something simply highlight the variable and click “>”. It will put the variable back in the master stored values box on the right.

Click Save and completely go out of the program.

Go back to the CDEMS opening screen > select a patient and do preview to see if the changes occurred.